

# MicroSENS CO<sub>2</sub> 180-HS 0-20%

## Incubator Sensor

CO<sub>2</sub>

Art. no. 7210.01-A.01



### NDIR gas sensor for CO<sub>2</sub> measurement in incubators, heat sterilizable at 180 °C.

The MicroSENS 180-HS 0-20% CO<sub>2</sub> gas sensor is an NDIR sensor for **measuring 5 vol.-% CO<sub>2</sub> in cell incubators to control ideal cell and tissue growth**. The sensor head is **placed in the incubation chamber** to measure gas concentration in the cell culture environment. Decontaminate the CO<sub>2</sub> diffusion sensor in an easy, effective and reliable way by **heat sterilization at typ. 180 °C**. The sensor head is designed for sterilization cycles in the four-digit range. It complies with DIN EN ISO 20857:2013, offering the highest number of sterilization cycles on the market. The heat sterilization with built-in sensor head eliminates the risk of cross-contamination. The sensor is **temperature and pressure compensated**. It features humidity compensation for best measurement results. Thanks to the **digital and analog interfaces**, the sensor is easy to integrate. Thanks to our **15-point calibration matrix** (five gas concentrations, at three temperatures each), the NDIR CO<sub>2</sub> sensor has an **improved linearity** compared to most sensors on the market. Achieve more precise measurements not only in the usual range of 5 vol.-%, but across the entire measurement range.

### Applications

- Medical diagnostics - in-vitro diagnostics
- Tests in cellular processes
- Biotechnology
- Pharmaceutical industry and research

### Product benefits

- High reliability
- Long service life
- Easy decontamination through high temperature sterilization cycle
- Low maintenance and minimal follow-up costs thanks to high-quality components

### Features

- IR dual-beam technology
- Temperature and pressure compensated
- Heat sterilizable at typ. 180 °C (maximum up to 190 °C)
- Humidity compensation
- Digital and analog interface
- Application-specific core components and electronic modules manufactured in-house
- Sensor head made of resistant and durable metal



## PRODUCT DATA SHEET

Gas Sensors / MicroSENS 180-HS 0-20% CO<sub>2</sub> gas sensor



## Technical data

### Measurement performance

Technical parameter	Value	Unit
Measuring gas	CO <sub>2</sub>	
Measuring range	0 ... 20	Vol.-%
Start-up time	< 15	s
Warm up time	< 5	min
IP protection class	IP20	
Repeatability	0.06	% FS
Accuracy	± (2 % of reading + 0.2 Vol.-%)	
Response time t90 with filter	< 45	s
Digital resolution	0.001	Vol.-%
Temperature influence	0.125	% FS / 10 K
Pressure influence	0.004	% FS / 10 hPa
Long-term stability	≤ ± 0.2 Vol.-% at 5 Vol.-% / year	

### Inputs and outputs

Operating voltage analog connection	15 ... 24	V DC
Operating voltage digital connection	12 ... 24	V DC
Nominal power consumption	< 2	W
Digital output signal	RS232	



## PRODUCT DATA SHEET

Gas Sensors / MicroSENS 180-HS 0-20% CO<sub>2</sub> gas sensor



### Operating environment

Technical parameter	Value	Unit
Operating temperature	0 ... 60	°C
Storage temperature	-20 ... 85	°C
Max. sterilization temperature <sup>1</sup>	190	°C
Humidity <sup>1</sup>	0 ... 95 % relative humidity (non condensing)	
Air pressure	600 ... 1200	hPa

<sup>1</sup> Sensor head only




## PRODUCT DATA SHEET

Gas Sensors / MicroSENS 180-HS 0-20% CO<sub>2</sub> gas sensor




## Typical operating characteristics

### Pin out




- Pin 1 – 24 V
- Pin 2 – RX
- Pin 3 – TX
- Pin 4 – GND
- Pin 5 – Current out

### Sensor overview




### Panel cut out



All geometrical dimensions in mm

### Combined panel cut out



## PRODUCT DATA SHEET

Gas Sensors / MicroSENS 180-HS 0-20% CO<sub>2</sub> gas sensor



## Disclaimer

All rights reserved. All information in this data sheet is based on latest knowledge, results of practical experience and tests carried out. Earlier specifications are hereby invalid. All specifications – technical included – are subject to change without notice. It is the customer's responsibility to ensure that the performance of the product is suitable for customer's specific application. No liability is accepted for indirect damage, in particular for the use or inability to use the product. Any liability we may have is limited to the value of the product itself.

